


# Introducing RaPET®—The Next Step. Stanbio's *Rapid Particle Enhanced Technology*.



**STANBIO**  
LABORATORY  
An EKF Diagnostics Company

## A complete line of highly sensitive and specific serology products.

Since 1960, Stanbio Laboratory has been Advancing the Science of Clinical Diagnostics. Our commitment to excellence ensures simply that "each and every test result is accurate each and every time." To this end we offer a complete line of highly sensitive and specific serology products.

Our line of RaPET® (Rapid Particle Enhanced Technology) latex products produces easy-to-read results in as little as 2 minutes. All procedures feature

"No Dilution" technology which eliminates the need to pre-dilute patient specimens. Use of color-coded dropper caps help to further reduce the chance of errors. All RaPET® assays utilize the same semi-quantitative testing procedure for the added convenience of using the same patient serial dilution for all RaPET® assays.

In addition, Stanbio offers a charcoal flocculation assay for the detection of Syphilis (RPR).

A Tri-level set of RPR/VDRL controls is available for VDRL users and other brands of RPR kits.

All Stanbio test kits come complete with test slide(s), pipettes / mixers, and controls to assure accurate test results. An easy to use worksheet for documenting semi-quantitative results is also included.

### The RaPET® Line

#### RaPET® ASO

- Latex agglutination slide test for the determination of Antistreptolysin-O in serum
- Obtain Qualitative or Semi-Quantitative results
- Positive/Negative Controls included
- Results in 2 minutes

50 Test  
100 Test

Ref. No. 1125-050  
Ref. No. 1125-100

#### RaPET® CRP

- Latex agglutination slide test for the determination of C-Reactive Protein in serum
- Obtain Qualitative or Semi-Quantitative results
- Positive/Negative Controls included
- Results in 2 minutes

50 Test  
100 Test

Ref. No. 1130-050  
Ref. No. 1130-100

#### RaPET® RF

- Latex agglutination slide test for the determination of Rheumatoid Factor in serum
- Obtain Qualitative or Semi-Quantitative results
- Positive/Negative Controls included
- Results in 2 minutes

50 Test  
100 Test

Ref. No. 1155-050  
Ref. No. 1155-100

#### RaPET® s-LE

- Latex agglutination slide test for the determination of Anti-DNP associated with Systemic Lupus Erythematosus in serum
- Obtain Qualitative or Semi-Quantitative results
- Positive/Negative Controls included
- Results in 3 minutes

50 Test

Ref. No. 1160-050



#### RPR Quicktest

- Carbon flocculation slide test for the determination of Reagin Antibodies associated with Syphilis in serum or plasma
- Obtain Qualitative or Quantitative results
- Controls Included—Reactive, Weakly Reactive and Non-Reactive
- Complete test kit includes pipettes / stirrers, dispensing needle, dispensing bottle and 18mm circle test cards
- Results in 8 minutes

100 Test  
500 Test

Ref. No. 1170-100  
Ref. No. 1170-500

#### Syphilis Tri-level Controls

- For use as quality control material in RPR/VDRL test procedures
- Liquid, ready-to-use controls
- Contains—Reactive, Weakly Reactive and Non-Reactive levels

3 x 5 mL

Ref. No. 1175-305



DISTRIBUTED BY:



*The World's Source For Clinical Diagnostics*

General & Special Chemistry / Near Patient Testing / Pregnancy / Serology  
Hematology / Occult Blood / Urinalysis / Rapid Tests  
Laboratory Detergent / Clinical Instrumentation / Pipettors

To order, or for more information,  
call: US Toll-Free: 866.STANBIO  
(782.6246)

Tel: 830.249.0772  
Fax: 830.249.0851

Stanbio Laboratory  
1261 North Main Street  
Boerne, Texas 78006 USA

stanbio@stanbio.com  
www.stanbio.com